

# ELISA LAB 1

Yellow pages and answers to post-lab questions are due May 22.

PROTOCOL I TRACKING DISEASE OUTBREAKS

## Laboratory Quick Guide


### Student Workstation Checklist

One workstation serves 4 students.

Item (Label)	Contents	Number	(✓)
Yellow tubes	Student test samples (0.75 ml)	4 (1 per student)	<input type="checkbox"/>
Violet tube (+)	Positive control (0.5 ml)	1	<input type="checkbox"/>
Blue tube (-)	Negative control (0.5 ml)	1	<input type="checkbox"/>
Green tube (PA)	Primary antibody (1.5 ml)	1	<input type="checkbox"/>
Orange tube (SA)	Secondary antibody (1.5 ml)	1	<input type="checkbox"/>
Brown tube (SUB)	Enzyme substrate (1.5 ml)	1	<input type="checkbox"/>
12-well microplate strips		2	<input type="checkbox"/>
50 $\mu$ l fixed-volume micropipet or 20–200 $\mu$ l adjustable micropipet		1	<input type="checkbox"/>
Yellow tips		10–20	<input type="checkbox"/>
Disposable plastic transfer pipets		5	<input type="checkbox"/>
70–80 ml wash buffer in beaker	Phosphate buffered saline with 0.05% Tween 20	1	<input type="checkbox"/>
Large stack of paper towels		2	<input type="checkbox"/>
Black marking pen		1	<input type="checkbox"/>

### ELISA for Tracking Disease Outbreaks

1. Label a yellow tube and a plastic transfer pipet with your initials.
2. Use the pipet to transfer all your "bodily fluid" sample into the tube of another student. Gently mix the samples, then take back half of the shared sample (750  $\mu$ l) to your own tube. Write down the name of the student next to "Sharing Partner #1".
3. When instructed to do so, repeat the sharing protocol two more times. Discard this transfer pipet after this step.



*Optional stopping point: Samples may be stored at 4°C overnight.*

It is ok if you write the person's # instead of name


Sharing Partner #1 \_\_\_\_\_

Sharing Partner #2 \_\_\_\_\_

Sharing Partner #3 \_\_\_\_\_

} write in lab notebook

4. Label your 12-well strip. On each strip label the first 3 wells with a "+" for the positive controls and the next 3 wells with a "-" for the negative controls. Label the remaining wells with your and your lab partner's initials (3 wells each).  
or yellow tube #




indicate in lab notebook who is Student A and who is student B

Use a sharpie to label




**PROTOCOL 1  
TRACKING DISEASE OUTBREAKS**

violet  
blue  
Your samples

5. Use a fresh pipet tip to transfer 50  $\mu$ l of the positive control (+) into the three "+" wells.
6. Use a fresh pipet tip to transfer 50  $\mu$ l of the negative control (-) into the three "-" wells.
7. Transfer 50  $\mu$ l of each of your team's samples from step 3 into the appropriately initialed three wells, using a fresh pipet tip for each sample.
8. Wait 5 minutes while all the proteins in the samples bind to the plastic wells.
9. WASH:



- a. Tip the microplate strip upside down onto the paper towels, and gently tap the strip a few times upside down. Make sure to avoid samples splashing back into wells.
- b. Discard the top paper towel.
- c. Use a fresh transfer pipet to fill each well with wash buffer, taking care not to spill over into wells. Note: the same transfer pipet is used for all washing steps.
- d. Tip the microplate strip upside down onto the paper towels and tap.
- e. Discard the top 2-3 paper towels.




wash buffer - ok to use straight out of bottle, don't need to put in beaker

- e. Discard the top 2-3 paper towels.  
(or replace them with fresh ones)
10. Repeat wash step 9.

**WASH**

11. Use a fresh pipet tip to transfer 50  $\mu$ l of primary antibody (PA) into all 12 wells of the microplate strip.
12. Wait 5 minutes for the antibodies to bind to their targets.
13. Wash the unbound primary antibody out of the wells by repeating all of wash step 9 **two** times.




**WASH 2x**

We may put paper towels atop of a sponge so we don't waste too many paper towels.

write in your lab notebook what you are doing as you do it

14. Use a fresh pipet tip to transfer 50  $\mu$ l of secondary antibody (SA) into all 12 wells of the microplate strip.




15. Wait 5 minutes for the antibodies to bind to their targets.


16. Wash the unbound secondary antibody out of the wells by repeating wash step 9 three times.

**WASH 3x**

17. Use a fresh pipet tip to transfer 50  $\mu$ l of enzyme substrate (SUB) into all 12 wells of the microplate strip.



18. Wait 5 minutes. Observe and record the results.



Record your results in your lab notebook.


## Student Manual

### Introduction

You are about to perform an experiment in which you will share simulated “body fluids” with your classmates. After sharing, you will perform an enzyme-linked immunosorbent assay or ELISA to determine if you have been exposed to a contagious “disease”. The ELISA uses antibodies to detect the presence of a disease agent, (for example, viruses, bacteria, or parasites) in your blood or other body fluid. You will then track the disease back to its source.

When you are exposed to a disease agent, your body mounts an immune response. Molecules that cause your body to mount an immune response are called antigens, and may include components of infectious agents like bacteria, viruses, and fungi. Within days, millions of antibodies — proteins that recognize the antigen and bind very tightly to it — are circulating in your bloodstream. Like magic bullets, antibodies seek out and attach themselves to their target antigens, flagging the invaders for destruction by other cells of the immune system.

Over 100 years ago, biologists found that animals’ immune systems respond to invasion by “foreign entities”, or antigens. Today, antibodies have become vital scientific tools, used in biotechnology research and to diagnose and treat disease. The number of different antibodies circulating in the blood has been estimated to be between  $10^6$  and  $10^{11}$ , so there is usually an antibody ready to deal with any antigen. In fact, antibodies make up to 15% of your total blood serum protein. Antibodies are very specific; each antibody recognizes only a single antigen.



A) Structure of IgG bound to the HIV capsid protein p24 as determined by X-ray crystallography (Harris et al. 1998, Momany et al. 1996). These structures can be downloaded from the Protein Data Bank ([www.pdb.uimg.br](http://www.pdb.uimg.br), (Berman et al. 2000) using the PDB identification codes 1IGY and 1AFV and manipulated using free online software such as Rasmol and Protein Explorer. B) A commonly used representation of an antibody bound to an antigen.

### How Are Antibodies Made?


Scientists have learned to use the immune response of animals to make antibodies that can be used as tools to detect and diagnose diseases. The study of the immune system is called “immunology”. Animals such as chickens, goats, rabbits, and sheep can be injected with an antigen and, after a period of time, their serum will contain antibodies that specifically recognize that antigen. If the antigen was a disease agent, the antibodies can be used to develop diagnostic tests for the disease. In an immunoassay, the antibodies used to recognize antigens like disease agents are called primary antibodies; primary antibodies confer specificity to the assay.

Other kinds of antibody tools, called secondary antibodies, are made in the same way. In an immunoassay, secondary antibodies recognize and bind to the primary antibodies, which are antibodies from another species. Secondary antibodies are prepared by injecting antibodies made in one species into another species. It turns out that antibodies from different species are different enough from each other that they will be recognized as foreign proteins and provoke an immune response. For example, to make a secondary antibody that will recognize a human primary antibody, human antibodies can be injected into an animal like a rabbit. After the rabbit mounts an immune response, the rabbit serum will contain antibodies that recognize and bind to human antibodies. The secondary antibodies used in this experiment are conjugated to the enzyme horseradish peroxidase (HRP) which produces a blue color in the presence of its substrate, TMB. These antibody and enzyme tools are the basis for the ELISA.

**Where Is ELISA Used in the Real World?**

With its rapid test results, the ELISA has had a major impact on many aspects of medicine and agriculture. ELISA is used for such diverse purposes as pregnancy tests, disease detection in people, animals, and plants, detecting illegal drug use, testing indoor air quality, and determining if food is labeled accurately. For new and emerging diseases like severe acute respiratory syndrome (SARS), one of the highest priorities of the US Centers for Disease Control (CDC) and the World Health Organization (WHO) has been to develop an ELISA that can quickly and easily verify whether patients have been exposed to the virus.

Some tests give positive or negative results in a matter of minutes. For example, home pregnancy dipstick tests are based on very similar principles to ELISA. They detect levels of human chorionic gonadotropin (hCG), a hormone that appears in the blood and urine of pregnant women within days of fertilization. The wick area of the dipstick is coated with anti-hCG antibody labeled with a pink compound (step 1). When the strip is dipped in urine, if hCG is present it will bind to the pink antibody, and the pink hCG-antibody complex will migrate up the strip via capillary action (step 2). When the pink complex reaches the first test zone, a narrow strip containing an unlabeled fixed anti-hCG antibody, the complex will bind and concentrate there, making a pink stripe (step 3). The dipsticks have a built-in control zone containing an unlabeled fixed secondary antibody that binds unbound pink complex (present in both positive and negative results) in the second stripe (step 4). Thus, every valid test will give a second pink stripe, but only a positive pregnancy test will give two pink stripes.



### **Why Do We Need Controls?**

Positive and negative controls are critical to any diagnostic test. Control samples are necessary to be sure your ELISA is working correctly. A positive control is a sample known to be positive for the disease agent, and a negative control is a sample that does not contain the disease agent.

### **Your Task Today**

You will be provided the tools and an experimental protocol to perform an ELISA. You will be given a simulated "body fluid" sample that you will share with your classmates. One or two of the samples in the class have been "infected". You will also be provided with positive and negative control samples. Then you and your fellow students will assay your samples for the presence of the "disease agent" to track the spread of the disease through your class population.

Now let's put this all together.

ELISA

~~ELISA~~ LAB 1

Answer on your own paper  
or in your lab notebook.  
due May 22 (yellow  
pages  
or answer  
answers on  
own paper.)

### Post-Lab Focus Questions

1. The samples that you added to the microplate strip contain many proteins and may or may not contain the disease antigen. What happened to the proteins in the plastic well if the sample contained the antigen? What if it did not contain the antigen?
2. Why did you need to wash the wells after every step?
3. When you added primary antibody to the wells, what happened if your sample contained the antigen? What if it did not contain the antigen?
4. When you added secondary antibody to the wells, what happened if your sample contained the antigen? What if it did not contain the antigen?
5. If the sample gave a negative result for the disease-causing agent, does this mean that you do not have the disease? What reasons could there be for a negative result when you actually do have the disease?
6. Why did you assay your samples in triplicate?
7. What antibody-based tests can you buy at your local pharmacy?
8. If you tested positive for disease exposure, did you have direct contact with one of the original infected students? If not, what conclusions can you reach about transmissibility of disease in a population?